

MiniLab – Net Force

1. An object of mass m is lifted by a person and accelerates upward amount a . Derive an expression for the force F applied by the person in terms of m , a , and any appropriate constants.
2. Test your result by using the Force and Acceleration sensor. Connect to Graphical Analysis – either by Bluetooth or USB connection. (see screenshots). Use ONLY sensor channels Force and x-axis Acceleration.
3. Once the sensor is connected, adjust the display to show two meters – notice the force and acceleration readouts. If held at rest on the table with the x-axis pointed upward, both should read zero (but probably will not). Click on the sensor control and Zero the acceleration. Also, still vertical and at rest, click to Zero the force reading as well (nothing should be touching the force “bumper”).
4. Once the sensors are zeroed, hold the sensor hanging by the force “bumper” (with a firm grip!). Click on the Collect button and then move the sensor smoothly up and down. Zoom in and inspect the graphs of Force and Acceleration vs. time – notice similarities? You should!
5. Create a single graph of force as a function of acceleration including line of best fit.
6. Attach the 50 gram mass to the bottom of the sensor using the thumbscrew. Collect new data as before – how does this affect the line of best fit?
7. Explain how the experiment supports or refutes the equation derived in question #1.
8. The slope equals what? The y-intercept equals what?



Vernier

Order Code: GDx-FOR

ID: 072051R8



Go Direct
Force and Acceleration

50 N
MAX

Untitled



Sensors



No Devices Connected

Connect to a wireless device below or connect via USB.

Discovered Wireless Devices

 PROXIMITY CONNECT

Filter Device List

e.g., 007 or TMP

 GDX-F0R 072051R8

Connect

DONE

y



Untitled

COLLECT



Sensors



Connected Devices

* GDX-FOR 072051R8 ⓘ

Disconnect

▼ SENSOR CHANNELS Force, X-axis acceleration

- Force
- X-axis acceleration
- Y-axis acceleration
- Z-axis acceleration
- X-axis gyro
- Y-axis gyro
- Z-axis gyro



Discovered Wireless Devices


PROXIMITY CONNECT

Searching for devices...

DONE

Force (N)


X-axis acceleration (m/s²)



Hold sensor at rest on table with the x-direction pointing vertically up and nothing touching the bumper.

Force and acceleration should both equal zero...

COLLECT



Hold sensor at rest on table with the x-direction pointing vertically up and nothing touching the bumper.

Force:
0.00 N

X-axis acceleration:
-0.05 m/s²

FORCE ⓘ
GDX-FOR 072051R8

Calibrate

Zero

Reverse

Units N ▼


X-AXIS ACCELERATION ⓘ
GDX-FOR 072051R8

Calibrate

Zero

Units m/s² ▼

COLLECT



Force:
0.88 N

X-axis
0.0



🔍

COLLECT



Graph Options



Title

Force vs. Acceleration

Appearance


 Points Lines Bars

x-axis range

Left to Right

Scaling Automatic

y-axis range

Bottom to Top 

view integral

Apply Curve Fit

Apply FFT

Apply Histogram


Add Annotation

Add Prediction

Edit Graph Options



Force vs. Acceleration



Force vs. Acceleration

