

Mini-Lab: $PV = NkT$

1. Connect GoDirect Pressure to Graphical Analysis on laptop. Move the plunger to 10 mL and then attach to the sensor by screwing it on snugly (not too tight!)
2. Under Data Collection Settings, choose Events with Entry with name Volume and units mL.
3. Click Collect and then Keep to record the current pressure. Enter the volume manually. Repeat the process to collect data for fixed volumes of 6, 8, 10, 12, 14, 16, 18, 20 mL.
4. Create a graph of Pressure vs. Volume. Determine the best fit. Use $PV = NkT$ to solve for the number of molecules in the syringe! Or, use 22.4 L to determine k .
5. Challenge analyze the same data with a linear graph.

Untitled

COLLECT



Data Collection Settings



Mode Event Based

Event Mode Events with Entry
 Selected Events


Event Name Volume

Units mL

 Average sensor reading over 10 seconds

CANCEL

DONE



Mode: Time Based Rate: 1 samples/s

Pressure: 99.4 kPa




Untitled

COLLECT

KEEP

(⌘) [Grid Icon] ...



Data Set 1		...	
Volume (mL)	...	Pressure (kPa)	...
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			



Mode: Events with Entry


Pressure: 99.3 kPa



STOP

KEEP

((a)) [] ...



Keep Point [X]

	Volume (mL)	Pressure (kPa)
1	<input type="text" value="10"/>	99.3

CANCEL **KEEP POINT**

Data Set 1 ...

Volume (mL) ...

Pressure (kPa) ...

1

2

3

4

5

6

7

8

9

10

